

KERN 327-11

Individual weight, OIML Class F1, 1 kg, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-441
Article number for DAkkS calibration	962-441
Article number for guaranteed plus-deviation	921-441

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	1 kg
Test weight - volume (cm ³)	125
Test weight - volume uncertainty (cm ³)	1,5625
Permissible tolerance according to OIML (+/-)	5 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	48×83 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	48×48×83 mm
Net weight	1 kg
Shipping method	Parcel service
Net weight approx.	1,00 kg
Gross weight approx.	1,00 kg
Shipping weight	1 kg