

KERN 327-13

Individual weight, OIML Class F1, 5 kg, Knob, stainless steel polished (OIML)



Services (optional)

Article number for conformity assessment	952-443
Article number for DAkkS calibration	962-443
Article number for guaranteed plus-deviation	921-443

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	625
Test weight - volume uncertainty (cm ³)	7,8125
Permissible tolerance according to OIML (+/-)	25 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	80×147 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	80×80×147 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg