

# KERN 327-131



Individual weight, OIML Class F1, 5 kg, Check weight, stainless steel polished (OIML)



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

## Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm <sup>3</sup> )	1250
Test weight - volume uncertainty (cm <sup>3</sup> )	15,625
Permissible tolerance according to OIML (+/-)	50 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

## Construction

Shape weight	Check weight
Dimensions (Ø×H)	136×132 mm
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	136×136×132 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

## Pictograms

### STANDARD

