

KERN 327-14



Individual weight, OIML Class F1, 10 kg, Knob, stainless steel polished (OIML)



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1250
Test weight - volume uncertainty (cm ³)	15,625
Permissible tolerance according to OIML (+/-)	50 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	100×185 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	100×100×185 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg

Services (optional)

Article number for conformity assessment	952-444
Article number for DAkkS calibration	962-444
Article number for guaranteed plus-deviation	921-444

Pictograms

STANDARD

