

KERN 327-151



Individual weight, OIML Class F1, 20 kg, Check weight, stainless steel polished (OIML)



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	20 kg
Test weight - volume (cm ³)	2500
Test weight - volume uncertainty (cm ³)	31,25
Permissible tolerance according to OIML (+/-)	100 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Check weight
Dimensions (Ø×H)	136×217 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	137×137×217 mm
Net weight	20 kg
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	20 kg
Shipping weight	20 kg

Services (optional)

Article number for DAkkS calibration	962-445
--------------------------------------	---------

Pictograms

STANDARD

