

KERN 327-16



Individual weight, OIML Class F1, 50 kg, Knob, stainless steel polished (OIML)



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F1
Product family	327

Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm ³)	6250
Test weight - volume uncertainty (cm ³)	78,125
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	8730
Weight OIML class	F1

Construction

Shape weight	Knob
Dimensions (Ø×H)	178×298 mm
Material	stainless steel polished (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	178×178×298 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Services (optional)

Article number for conformity assessment	952-446
Article number for DAkkS calibration	962-446
Article number for guaranteed plus-deviation	921-446

Pictograms

STANDARD

