

# KERN 328-22

**KERN**

Set of milligram weights, F1, 1 mg - 500 mg, Platelet, stainless steel polished (OIML), in plastic case



Article number for guaranteed plus-deviation 921-450

Pictograms

STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F1
Product family	328

## Measuring System

Test weight - nominal value	1 mg - 500 mg
Weight OIML class	F1

## Construction

Shape weight	Platelet
Dimension (W×D×H)	136×53×22 mm
Number of weights	13
Case of weight set	in plastic case
Material	stainless steel polished (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	130×60×22 mm
Net weight	0,1 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

## Services (optional)

Article number for conformity assessment	952-450
Article number for DAkkS calibration	962-450