

KERN 328-222-410

KERN

Plastic box for weight sets (Cylindrical/compact, 1 mg - 5 g)



Category	
Brand	KERN
Product category	Cases & Boxes
Product group	Plastic box

Measuring System	
Case suitable for OIML classes	E1 - M2
Case suitable for nominal weight values	1 mg - 5 g
Suitable weight shape for case	Cylindrical/compact

Construction	
Dimension (W×D×H)	130×60×22 mm
Material	plastic

Packing & Shipping	
Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	130×60×22 mm
Net weight	0,09 kg
Shipping method	Parcel service
Net weight approx.	0,10 kg
Gross weight approx.	0,10 kg
Shipping weight	0,1 kg

Pictograms	
------------	--