

# KERN 333-036

# KERN

Weight set, F2, 1 mg - 100 g, Knob, stainless steelfine turned (OIML), in aluminium case



## Pictograms

### STANDARD



### Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

### Measuring System

Test weight - nominal value	1 mg - 100 g
Weight OIML class	F2

### Construction

Shape weight	Knob
Dimension (W×D×H)	165×125×75 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

### Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×125×75 mm
Net weight	0,65 kg
Shipping method	Parcel service
Net weight approx.	0,65 kg
Gross weight approx.	0,90 kg
Shipping weight	0,9 kg

### Services (optional)

Article number for DAkkS calibration	962-402
Article number for guaranteed plus-deviation	921-402