

# KERN 333-04

**KERN**

Weight set, F2, 1 mg - 200 g, Knob, stainless steelfine turned (OIML), in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

## Measuring System

Test weight - nominal value	1 mg - 200 g
Weight OIML class	F2

## Construction

Shape weight	Knob
Dimension (W×D×H)	160×98×58 mm
Number of weights	23
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	160×98×58 mm
Net weight	1,17 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	1,2 kg
Shipping weight	1,17 kg

## Services (optional)

Article number for DAkkS calibration	962-403
Article number for guaranteed plus-deviation	921-403