

KERN 333-046

KERN

Weight set, F2, 1 mg - 200 g, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

Measuring System

Test weight - nominal value	1 mg - 200 g
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	190×130×80 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

Services (optional)

Article number for DAkkS calibration	962-403
Article number for guaranteed plus-deviation	921-403