

KERN 333-06

KERN

Weight set, F2, 1 mg - 1 kg, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	230×180×112 mm
Number of weights	25
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	230×180×112 mm
Net weight	4,21 kg
Shipping method	Parcel service
Net weight approx.	4,4 kg
Gross weight approx.	4,4 kg
Shipping weight	4,21 kg

Services (optional)

Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405