

# KERN 333-064

**KERN**

Weight set, F2, 1 mg - 1 kg, Knob, stainless steelfine turned (OIML), in plastic case



## Pictograms

### STANDARD



## Category

|                  |            |
|------------------|------------|
| Brand            | KERN       |
| Product category | Weight set |
| Product group    | OIML F2    |
| Product family   | 333        |

## Measuring System

|                             |             |
|-----------------------------|-------------|
| Test weight - nominal value | 1 mg - 1 kg |
| Weight OIML class           | F2          |

## Construction

|                    |                                   |
|--------------------|-----------------------------------|
| Shape weight       | Knob                              |
| Dimension (W×D×H)  | 350×320×180 mm                    |
| Number of weights  | 25                                |
| Case of weight set | in plastic case                   |
| Material           | stainless steelfine turned (OIML) |

## Packing & Shipping

|                              |                |
|------------------------------|----------------|
| Readability force [d] (N)    | 1 d            |
| Dimensions packaging (W×D×H) | 350×320×180 mm |
| Net weight                   | 3,7 kg         |
| Shipping method              | Parcel service |
| Net weight approx.           | 3,8 kg         |
| Gross weight approx.         | 3,8 kg         |
| Shipping weight              | 4 kg           |

## Services (optional)

|  |         |
|--|---------|
| Article number for DAkkS calibration         | 962-405 |
| Article number for guaranteed plus-deviation | 921-405 |