

KERN 333-066

KERN

Weight set, F2, 1 mg - 1 kg, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

Measuring System

Test weight - nominal value	1 mg - 1 kg
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	2,7 kg
Shipping method	Parcel service
Net weight approx.	2,8 kg
Gross weight approx.	3,6 kg
Shipping weight	4 kg

Services (optional)

Article number for DAkkS calibration	962-405
Article number for guaranteed plus-deviation	921-405