

KERN 333-084

KERN

Weight set, F2, 1 mg - 5 kg, Knob, stainless steelfine turned (OIML), in plastic case



Article number for guaranteed
plus-deviation

921-407

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

Measuring System

Test weight - nominal value	1 mg - 5 kg
OIML-compliant	✓
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	28
Case of weight set	in plastic case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	12,69 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	13 kg
Shipping weight	12,69 kg

Services (optional)

Article number for DAKkS calibration	962-407
---	---------