

KERN 333-086

KERN

Weight set, F2, 1 mg - 5 kg, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	333

Measuring System

Test weight - nominal value	1 mg - 5 kg
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	460×260×200 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Net weight	12,6 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	18 kg
Shipping weight	18 kg

Services (optional)

Article number for DAkkS calibration	962-407
Article number for guaranteed plus-deviation	921-407