

KERN 334-02

KERN

Weight set, F2, 1 g - 50 g, Knob, stainless steelfine turned (OIML), in wooden case



Article number for guaranteed
plus-deviation

921-415

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

Measuring System

Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	129×86×40 mm
Number of weights	8
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	129×86×40 mm
Net weight	0,26 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,26 kg

Services (optional)

Article number for DAKkS calibration	962-415
---	---------