

KERN 334-04

KERN

Weight set, F2, 1 g - 200 g, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	160×98×58 mm
Number of weights	11
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	160×98×58 mm
Net weight	1,04 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	1,2 kg
Shipping weight	1,04 kg

Services (optional)

Article number for DAkkS calibration	962-417
Article number for guaranteed plus-deviation	921-417