

# KERN 334-064

# KERN

Weight set, F2, 1 g - 1 kg, Knob, stainless steelfine turned (OIML), in plastic case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

## Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	F2

## Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	13
Case of weight set	in plastic case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	3,7 kg
Shipping method	Parcel service
Net weight approx.	3,8 kg
Gross weight approx.	3,8 kg
Shipping weight	4 kg

## Services (optional)

Article number for DAkkS calibration	962-419
Article number for guaranteed plus-deviation	921-419