

KERN 334-07

KERN

Weight set, F2, 1 g - 2 kg, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML F2 |
| Product family | 334 |

Measuring System

| | |
|-----------------------------|------------|
| Test weight - nominal value | 1 g - 2 kg |
| Weight OIML class | F2 |

Construction

| | |
|--------------------|-----------------------------------|
| Shape weight | Knob |
| Dimension (WxDxH) | 305x180x140 mm |
| Number of weights | 15 |
| Case of weight set | in wooden case |
| Material | stainless steelfine turned (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (WxDxH) | 305x180x140 mm |
| Net weight | 9,01 kg |
| Shipping method | Parcel service |
| Net weight approx. | 10 kg |
| Gross weight approx. | 10 kg |
| Shipping weight | 9,01 kg |

Services (optional)

| | |
|--|---------|
| Article number for DAkkS calibration | 962-420 |
| Article number for guaranteed plus-deviation | 921-420 |