

KERN 334-074

KERN

Weight set, F2, 1 g - 2 kg, Knob, stainless steelfine turned (OIML), in plastic case



Article number for guaranteed
plus-deviation

921-420

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

Measuring System

Test weight - nominal value	1 g - 2 kg
OIML-compliant	✓
Weight OIML class	F2

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	15
Case of weight set	in plastic case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	7,69 kg
Shipping method	Parcel service
Net weight approx.	8 kg
Gross weight approx.	8 kg
Shipping weight	7,69 kg

Services (optional)

Article number for DAkkS calibration	962-420
---	---------