

# KERN 334-08

**KERN**

Weight set, F2, 1 g - 5 kg, Knob, stainless steelfine turned (OIML), in wooden case



Article number for guaranteed  
plus-deviation

921-421

## Pictograms

## STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

## Measuring System

Test weight - nominal value	1 g - 5 kg
OIML-compliant	✓
Weight OIML class	F2

## Construction

Shape weight	Knob
Dimension (W×D×H)	430×190×200 mm
Number of weights	16
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	430×190×200 mm
Net weight	17,81 kg
Shipping method	Parcel service
Net weight approx.	18 kg
Gross weight approx.	18 kg
Shipping weight	17,81 kg

## Services (optional)

Article number for DAkkS calibration	962-421
---	---------