

# KERN 334-086

**KERN**

Weight set, F2, 1 g - 5 kg, Knob, stainless steelfine turned (OIML), in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML F2
Product family	334

## Measuring System

Test weight - nominal value	1 g - 5 kg
Weight OIML class	F2

## Construction

Shape weight	Knob
Dimension (W×D×H)	460×260×200 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	460×260×200 mm
Net weight	12,6 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	18 kg
Shipping weight	18 kg

## Services (optional)

Article number for DAkkS calibration	962-421
Article number for guaranteed plus-deviation	921-421