

# KERN 337-1 10-200

**KERN**

Wooden box for individual weights (Cylindrical/compact, 1 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	1 kg
Suitable weight shape for case	Cylindrical/compact

## Construction

Dimension (W×D×H)	73×73×93 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	73×73×93 mm
Net weight	0,22 kg
Shipping method	Parcel service
Net weight approx.	0,25 kg
Gross weight approx.	0,25 kg
Shipping weight	0,235 kg

## Pictograms

STANDARD

