

KERN 337-130-200

KERN

Wooden box for individual weights (Cylindrical/compact, 5 kg)



Category

| | |
|------------------|---------------|
| Brand | KERN |
| Product category | Cases & Boxes |
| Product group | Wooden box |

Measuring System

| | |
|---|---------------------|
| Case suitable for OIML classes | F2 + M1 |
| Case suitable for nominal weight values | 5 kg |
| Suitable weight shape for case | Cylindrical/compact |

Construction

| | |
|-------------------|----------------|
| Dimension (W×D×H) | 106×106×170 mm |
| Material | beech |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 106×106×170 mm |
| Net weight | 0,52 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,55 kg |
| Gross weight approx. | 0,55 kg |
| Shipping weight | 0,52 kg |

Pictograms

STANDARD

