

# KERN 337-141-100

**KERN**

Wooden box for individual weights (Check weight, 10 kg)



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

## Measuring System

Case suitable for OIML classes	F1
Case suitable for nominal weight values	10 kg
Suitable weight shape for case	Check weight

## Construction

Dimension (W×D×H)	205×205×190 mm
Material	beech

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	205×205×190 mm
Net weight	3,1 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,2 kg
Shipping weight	3,1 kg

## Pictograms

STANDARD

