

KERN 337-141

KERN

Individual weight, OIML Class F2, 10 kg, Check weight, stainless steelfine turned (OIML)



Services (optional)

Article number for DAkkS calibration	962-444
Article number for guaranteed plus-deviation	921-444

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Test weight - nominal value	10 kg
Test weight - volume (cm ³)	1250
Test weight - volume uncertainty (cm ³)	15,625
Permissible tolerance according to OIML (+/-)	160 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	10700
Weight OIML class	F2

Construction

Shape weight	Check weight
Dimensions (Ø×H)	136×132 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	136×136×132 mm
Net weight	10 kg
Shipping method	Parcel service
Net weight approx.	10 kg
Gross weight approx.	10 kg
Shipping weight	10 kg