

KERN 337-150-200

KERN

Wooden box for individual weights (Cylindrical/compact, 20 kg)



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Wooden box

Measuring System

Case suitable for OIML classes	F2 + M1
Case suitable for nominal weight values	20 kg
Suitable weight shape for case	Cylindrical/compact

Construction

Dimension (W×D×H)	185×185×300 mm
Material	beech

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	185×185×300 mm
Net weight	5,2 kg
Shipping method	Parcel service
Net weight approx.	6 kg
Gross weight approx.	6 kg
Shipping weight	5,2 kg

Pictograms

STANDARD

