

KERN 337-161



Individual weight, OIML Class F2, 50 kg, Check weight, stainless steelfine turned (OIML)



Category

Brand	KERN
Product category	Individual weight
Product group	OIML F2
Product family	337

Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm ³)	6250
Test weight - volume uncertainty (cm ³)	78,125
Permissible tolerance according to OIML (+/-)	800 mg
Test weight - Density [Min] (OIML)	7390
Test weight - Density [Max] (OIML)	10700
Weight OIML class	F2

Construction

Shape weight	Check weight
Dimensions (Ø×H)	198×249,5 mm
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	198×198×249 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

Services (optional)

Article number for DAkkS calibration	962-446
Article number for guaranteed plus-deviation	921-446

Pictograms

STANDARD

