

KERN 343-03

KERN

Weight set, M1, 1 mg - 100 g, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

| | |
|------------------|------------|
| Brand | KERN |
| Product category | Weight set |
| Product group | OIML M1 |
| Product family | 343 |

Measuring System

| | |
|-----------------------------|--------------|
| Test weight - nominal value | 1 mg - 100 g |
| Weight OIML class | M1 |

Construction

| | |
|--------------------|-----------------------------------|
| Shape weight | Knob |
| Dimension (W×D×H) | 151×91×49 mm |
| Number of weights | 21 |
| Case of weight set | in wooden case |
| Material | stainless steelfine turned (OIML) |

Packing & Shipping

| | |
|------------------------------|----------------|
| Readability force [d] (N) | 1 d |
| Dimensions packaging (W×D×H) | 151×91×49 mm |
| Net weight | 0,63 kg |
| Shipping method | Parcel service |
| Net weight approx. | 0,65 kg |
| Gross weight approx. | 0,65 kg |
| Shipping weight | 0,63 kg |

Services (optional)

| | |
|--|---------|
| Article number for conformity assessment | 952-602 |
| Article number for DAkkS calibration | 962-602 |