

KERN 343-04

KERN

Weight set, M1, 1 mg - 200 g, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 200 g
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	160×98×58 mm
Number of weights	23
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	160×98×58 mm
Net weight	1,17 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	1,2 kg
Shipping weight	1,17 kg

Services (optional)

Article number for conformity assessment	952-603
Article number for DAkkS calibration	962-603