

# KERN 343-05

# KERN

Weight set, M1, 1 mg - 500 g, Knob, stainless steelfine turned (OIML), in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

## Measuring System

Test weight - nominal value	1 mg - 500 g
Weight OIML class	M1

## Construction

Shape weight	Knob
Dimension (W×D×H)	168×98×78 mm
Number of weights	24
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	168×98×78 mm
Net weight	1,76 kg
Shipping method	Parcel service
Net weight approx.	1,8 kg
Gross weight approx.	1,8 kg
Shipping weight	1,76 kg

## Services (optional)

Article number for conformity assessment	952-604
Article number for DAKkS calibration	962-604