

KERN 343-054

KERN

Weight set, M1, 1 mg - 500 g, Knob, stainless steelfine turned (OIML), in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 500 g
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	250×210×75 mm
Number of weights	24
Case of weight set	in plastic case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	250×210×75 mm
Net weight	1,44 kg
Shipping method	Parcel service
Net weight approx.	1,6 kg
Gross weight approx.	1,6 kg
Shipping weight	1,44 kg

Services (optional)

Article number for conformity assessment	952-604
Article number for DAKkS calibration	962-604