

KERN 343-056

KERN

Weight set, M1, 1 mg - 500 g, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 500 g
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	210×140×98 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	210×140×98 mm
Net weight	1,61 kg
Shipping method	Parcel service
Net weight approx.	1,8 kg
Gross weight approx.	2,0 kg
Shipping weight	2 kg

Services (optional)

Article number for conformity assessment	952-604
Article number for DAKkS calibration	962-604