

KERN 343-08

KERN

Weight set, M1, 1 mg - 5 kg, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 5 kg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	430×190×200 mm
Number of weights	28
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	430×190×200 mm
Net weight	17,81 kg
Shipping method	Parcel service
Net weight approx.	18 kg
Gross weight approx.	18 kg
Shipping weight	17,81 kg

Services (optional)

Article number for conformity assessment	952-607
Article number for DAkkS calibration	962-607