

KERN 343-084

KERN

Weight set, M1, 1 mg - 5 kg, Knob, stainless steelfine turned (OIML), in plastic case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 5 kg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	350×320×180 mm
Number of weights	28
Case of weight set	in plastic case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	350×320×180 mm
Net weight	12,69 kg
Shipping method	Parcel service
Net weight approx.	13 kg
Gross weight approx.	13 kg
Shipping weight	12,69 kg

Services (optional)

Article number for conformity assessment	952-607
Article number for DAkkS calibration	962-607