

KERN 343-09

KERN

Weight set, M1, 1 mg - 10 kg, Knob, stainless steelfine turned (OIML), in wooden case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	540×190×245 mm
Number of weights	29
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	540×190×245 mm
Net weight	30,61 kg
Shipping method	Freight
Net weight approx.	32 kg
Gross weight approx.	34 kg
Shipping weight	32,5 kg

Services (optional)

Article number for conformity assessment	952-608
Article number for DAkkS calibration	962-608