

KERN 343-096

KERN

Weight set, M1, 1 mg - 10 kg, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	343

Measuring System

Test weight - nominal value	1 mg - 10 kg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	610×260×235 mm
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	610×260×235 mm
Net weight	23,91 kg
Shipping method	Parcel service
Net weight approx.	24 kg
Gross weight approx.	28 kg
Shipping weight	28 kg

Services (optional)

Article number for conformity assessment	952-608
Article number for DAKkS calibration	962-608