

KERN 344-02

Weight set, M1, 1 g - 50 g, Knob, stainless steelfine turned (OIML), in wooden case



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

Measuring System

Test weight - nominal value	1 g - 50 g
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	129×86×40 mm
Number of weights	8
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	129×86×40 mm
Net weight	0,26 kg
Shipping method	Parcel service
Net weight approx.	0,30 kg
Gross weight approx.	0,30 kg
Shipping weight	0,26 kg

Services (optional)

Article number for conformity assessment	952-615
Article number for DAkkS calibration	962-615

Pictograms

STANDARD

