

KERN 344-026

KERN

Weight set, M1, 1 g - 50 g, Knob, stainless steelfine turned (OIML), in aluminium case



Article number for DAkkS
calibration

962-615

Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

Measuring System

Test weight - nominal value	1 g - 50 g
OIML-compliant	✓
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	165×120×75 mm
Number of weights	8
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Net weight	0,55 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

Services (optional)

Article number for conformity assessment	952-615
---	---------