

# KERN 344-026

# KERN

Weight set, M1, 1 g - 50 g, Knob, stainless steelfine turned (OIML), in aluminium case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

## Measuring System

Test weight - nominal value	1 g - 50 g
Weight OIML class	M1

## Construction

Shape weight	Knob
Dimension (W×D×H)	165×120×75 mm
Number of weights	8
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	165×120×75 mm
Net weight	0,55 kg
Shipping method	Parcel service
Net weight approx.	0,55 kg
Gross weight approx.	0,80 kg
Shipping weight	0,8 kg

## Services (optional)

Article number for conformity assessment	952-615
Article number for DAkkS calibration	962-615