

# KERN 344-04

# KERN

Weight set, M1, 1 g - 200 g, Knob, stainless steelfine turned (OIML), in wooden case



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

## Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	M1

## Construction

Shape weight	Knob
Dimension (W×D×H)	160×98×58 mm
Number of weights	11
Case of weight set	in wooden case
Material	stainless steelfine turned (OIML)

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	160×98×58 mm
Net weight	1,04 kg
Shipping method	Parcel service
Net weight approx.	1,2 kg
Gross weight approx.	1,2 kg
Shipping weight	1,04 kg

## Services (optional)

Article number for conformity assessment	952-617
Article number for DAkkS calibration	962-617