

KERN 344-046

KERN

Weight set, M1, 1 g - 200 g, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

Measuring System

Test weight - nominal value	1 g - 200 g
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (W×D×H)	190×130×80 mm
Number of weights	11
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	190×130×80 mm
Net weight	0,9 kg
Shipping method	Parcel service
Net weight approx.	0,90 kg
Gross weight approx.	1,2 kg
Shipping weight	1,1 kg

Services (optional)

Article number for conformity assessment	952-617
Article number for DAKkS calibration	962-617