

KERN 344-066

KERN

Weight set, M1, 1 g - 1 kg, Knob, stainless steelfine turned (OIML), in aluminium case



Pictograms

STANDARD



Category

Brand	KERN
Product category	Weight set
Product group	OIML M1
Product family	344

Measuring System

Test weight - nominal value	1 g - 1 kg
Weight OIML class	M1

Construction

Shape weight	Knob
Dimension (WxDxH)	230x140x110 mm
Number of weights	13
Case of weight set	in aluminium case
Material	stainless steelfine turned (OIML)

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (WxDxH)	230x140x110 mm
Net weight	3,15 kg
Shipping method	Parcel service
Net weight approx.	3,2 kg
Gross weight approx.	3,6 kg
Shipping weight	3,5 kg

Services (optional)

Article number for conformity assessment	952-619
Article number for DAKkS calibration	962-619