


KERN 346-022-005

KERN

Weight container for rectangular weights, calibrated to OIML class M1



Pictograms

STANDARD



Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container

Measuring System

Test weight - nominal value	20 kg
Permissible tolerance according to OIML (+/-)	1 mg
Weight OIML class	M1

Construction

Shape weight	weight container
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Net weight	20 kg
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	20 kg
Shipping weight	50 kg

Services (optional)

Article number for DAkkS calibration	962-645
--------------------------------------	---------