


# KERN 346-052-010

**KERN**

Weight container for rectangular weights, calibrated to OIML class M1



## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Cases & Boxes
Product group	Weight container

## Measuring System

Test weight - nominal value	50 kg
Permissible tolerance according to OIML (+/-)	2,5 mg
OIML-compliant	✓
Weight OIML class	M1

## Construction

Shape weight	weight container
Material	stainless steel

## Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	1000×500×500 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg

## Services (optional)

Article number for DAkkS calibration	962-646
--------------------------------------	---------