

KERN 346-06

KERN

Individual weight, OIML Class M1, 5 kg, Block, stainless steel



Services (optional)

Article number for conformity assessment	952-643
Article number for DAkkS calibration	962-643

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	5 kg
Test weight - volume (cm ³)	694,4
Test weight - volume uncertainty (cm ³)	57,9
Permissible tolerance according to OIML (+/-)	250 mg
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Block
Dimensions weight (W×D×H)	150×68×96 mm
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	150×68×96 mm
Net weight	5 kg
Shipping method	Parcel service
Net weight approx.	5,0 kg
Gross weight approx.	5,0 kg
Shipping weight	5 kg