

# KERN 346-08

**KERN**

Individual weight, OIML Class M1, 20 kg, Block, stainless steel



## Services (optional)

Article number for conformity assessment	952-645
Article number for DAkkS calibration	962-645

## Pictograms

### STANDARD



## Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

## Measuring System

Test weight - nominal value	20 kg
Test weight - volume (cm <sup>3</sup> )	2778
Test weight - volume uncertainty (cm <sup>3</sup> )	231
Permissible tolerance according to OIML (+/-)	1 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

## Construction

Shape weight	Block
Dimensions weight (W×D×H)	233×115×147 mm
Material	stainless steel

## Packing & Shipping

Readability force [d] (N)	1 d
Dimensions packaging (W×D×H)	233×115×145 mm
Net weight	20 kg
Shipping method	Parcel service
Net weight approx.	20 kg
Gross weight approx.	20 kg
Shipping weight	20 kg