

KERN 346-09

KERN

Individual weight, OIML Class M1, 50 kg, Block, stainless steel



Services (optional)

Article number for conformity assessment	952-646
Article number for DAkkS calibration	962-646

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	50 kg
Test weight - volume (cm ³)	6944
Test weight - volume uncertainty (cm ³)	579
Permissible tolerance according to OIML (+/-)	2,5 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Block
Dimensions weight (W×D×H)	295×158×191 mm
Material	stainless steel

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	295×158×191 mm
Net weight	50 kg
Shipping method	Freight
Net weight approx.	50 kg
Gross weight approx.	50 kg
Shipping weight	50 kg