

KERN 346-82

Individual weight, OIML Class M1, 200 kg, Heavy Duty Weight, cast iron lacquered



Services (optional)

Article number for DAkkS calibration	962-692
--------------------------------------	---------

Pictograms

STANDARD



Category

Brand	KERN
Product category	Individual weight
Product group	OIML M1
Product family	346

Measuring System

Test weight - nominal value	200 kg
Test weight - volume (cm ³)	27778
Test weight - volume uncertainty (cm ³)	2315
Permissible tolerance according to OIML (+/-)	10 g
Test weight - Density [Min] (OIML)	3000
Weight OIML class	M1

Construction

Shape weight	Heavy Duty Weight
Dimensions weight (W×D×H)	465×340×291 mm
Material	cast iron lacquered

Packing & Shipping

Readability force [d] (N)	2 d
Dimensions packaging (W×D×H)	465×340×291 mm
Net weight	200 kg
Shipping method	Freight
Net weight approx.	200 kg
Gross weight approx.	200 kg
Shipping weight	200 kg